


Dental Diamond Burs Catalog

High Quality Diamond Dental Instruments

CE 0197

Contact Information

Email: dentallabshop@163.com

WhatsApp: +86 15514831900




Dental Diamond Bur

牙科金刚砂车针

Advantages

1. Strictly comply with ISO Standard ;
2. Sharper and More Efficient ;
3. Durable and Valuable ;
4. Full Range of shapes and grits.
5. Fastly develop new items.
6. Individual Design and OEM Service



001-008M

Item No.

Grit

- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine
- Yellow / Extra Fine

0.8 19.0 001(801) BR-49	1.0 19.0 001(801) BR-45	1.2 19.0 001(801) BR-46	1.4 19.0 001(801) BR-41	1.6 19.0 001(801) BR-40	1.8 19.0 001(801) BR-31	2.1 19.0 001(801) BR-28	2.3 19.0 001(801) BR-26	2.5 19.0 001(801) BR-25	2.9 19.0 001(801) BR-23
001-008M	001-010M	001-012M	001-014M	001-016M	001-018M	001-021M	001-023M	001-025M	001-029M

圆球型 Ball, Round /BR

6mm 1.2 22.0 697(801L) BR-L46	1.4 22.0 697(801L) BR-L41	1.6 22.0 697(801L) BR-L40	1.8 22.0 697(801L) BR-L31	1.2/2.3 19.0 002(802) BC-43	1.4/2.3 19.0 002(802) BC-42	1.6/3.5 19.0 002(802) BC-32	1.8/2.3 19.0 002(802) BC-31	1.2/10.0 21.0 494(802L) ENDO-012	1.4/10.0 21.0 494(802L) ENDO-014
697-012M	697-014M	697-016M	697-018M	002-012M	002-014M	002-016M	002-018M	494-012M	494-014M

圆球型 Round with Long Neck /BR

紧口圆球型 Round with Collor /BC

1.6/10.0 21.0 494(802L) ENDO-016	1.8/10.0 21.0 494(802L) ENDO-018	1.0(0.8)/1.0 19.0 010(805) SI-45	1.2(0.9)/1.2 19.0 010(805) SI-46	1.4(1.0)/1.6 19.0 010(805) SI-47	1.6(1.1)/1.8 19.0 010(805) SI-48	1.4(1.1)/4.0 19.0 225(807) SI-50	1.6(1.2)/4.0 19.0 225(807) ---	1.8(1.3)/5.0 19.0 225(807) SI-51
494-016M	494-018M	010-010M	010-012M	010-014M	010-016M	225-014M	225-016M	225-018M

倒锥型 Inverted Cone / SI

1.0/2.3 19.0 019(806) DI-41	1.4/3.2 19.0 019(806) DI-42	3.2/5.0 19.0 038(811) EX-11	3.5/7.0 20.0 039(811L) EX-12	4.5/0.7 19.0 044(818) WR-44	4.2/0.9 19.0 067(909) WR-11	4.2/1.8 19.0 068(909) WR-13
019-010M	019-014M	038-032M	039-035M	044-045M	067-042M	068-042M

紧口倒锥型 Inverted Cone with collar/DI

双圆锥型 Double cone/EX

轮型 Wheel/WR

1.0/4.0 19.0 109(835) SF-41	1.2/4.0 19.0 109(835) SF-31	1.4/4.0 19.0 109(835) SF-42	1.6/4.0 19.0 109(835) ---	1.2/6.0 21.0 110(836) SF-20	1.4/6.0 21.0 110(836) SF-21	1.6/6.0 21.0 110(836) SF-15	1.2/8.0 22.0 111(837) SF-11	1.4/8.0 22.0 111(837) SF-12	1.6/8.0 22.0 111(837) SF-13
109-010M	109-012M	109-014M	109-016M	110-012M	110-014M	110-016M	111-012M	111-014M	111-016M

平头圆柱型 Flat End Cylinder/SF

1.2/8.0 22.0 130(885) CP-11	1.4/8.0 22.0 130(885) CP-12	1.6/8.0 22.0 130(885) CP-13	1.2/10.0 24.0 131(886) ---	1.4/10.0 24.0 131(886) ---	1.6/10.0 24.0 131(886) ---	1.4(0.7)/10.0 24.0 173(848) TF-11
130-012M	130-014M	130-016M	131-012M	131-014M	131-016M	173-014M

尖头圆柱型 Ogival end Cylinder/CP

173-014M — 规格
平头锥型 — 型号

金刚石粒度

- 黑色/极粗粒度
- 绿色/粗粒度
- 蓝色/标准粒度
- 红色/较细粒度
- 黄色/极细粒度

	1.0/4.0 19.0 139(838) SR-41		1.2/4.0 19.0 139(838) SR-31		1.4/4.0 19.0 139(838) SR-42		1.4/6.0 21.0 140(880) SR-21		1.2/8.0 22.0 141(881) SR-11		1.4/8.0 22.0 141(881) SR-12		1.6/8.0 22.0 141(881) SR-13		1.2/0.2 21.0 150(839) EX-18
139-010M 139-012M 139-014M 140-014M 141-012M 141-014M 141-016M 150-012M															
圆头柱型 Round End, Cylinder/SR												肩台车针 End Cut			

	1.2/6.0 21.0 164(852) TC-26		1.4/8.0 22.0 165(858) TC-21		1.0/10.0 22.0 166(859) TC-27		1.2/10.0 22.0 166(859) TC-28		1.4/10.0 22.0 166(859) TC-10		1.6/10.0 22.0 166(859) TC-11		1.8/10.0 22.0 166(859) TC-16
164-012M 165-014M 166-010M 166-012M 166-014M 166-016M 166-018M													
尖头锥型 Pointed cone/TC													



	1.0(0.6)/3.0 19.0 169(845) TF-41		1.2(0.7)/4.0 19.0 170(845) TF-42		1.4(0.8)/4.0 19.0 170(845) TF-43		1.6(0.9)/4.0 19.0 170(845) TF-31		1.8(1.0)/5.0 20.0 170(845) TF-23		2.1(1.5)/5.0 20.0 170(845) TF-22		1.4(0.7)/6.5 21.0 171(846) TF-20		1.6(0.8)/6.5 21.0 171(846) TF-21		1.0(0.6)/8.0 22.0 172(847) TF-07		1.2(0.7)/8.0 22.0 172(847) TF-08
169-010M 170-012M 170-014M 170-016M 170-018M 170-021M 171-014M 171-016M 172-010M 172-012M																			
平头锥型 Flat End Cone/TF																			

	1.4(0.7)/8.0 22.0 172(847) TF-09		1.6(0.8)/8.0 22.0 172(847) TF-10		1.8(0.9)/8.0 22.0 172(847) TF-56		2.3(1.3)/9.0 22.0 172(847) TF-14		1.4(0.7)/10.0 24.0 173(848) TF-11		1.6(0.8)/10.0 22.0 173(848) TF-12		1.8(0.9)/10.0 22.0 173(848) TF-13
172-014M 172-016M 172-018M 172-023M 173-014M 173-016M 173-018M													

	1.0(0.6)/4.0 19.0 196(849) CR-41		1.2(0.7)/4.0 19.0 196(849) CR-31		1.4(0.8)/4.0 19.0 196(849) CR-42		1.6(0.9)/4.0 19.0 196(849) CR-12		2.0(1.0)/5.0 19.0 196(849) CR-11		1.2(0.7)/6.5 21.0 197(855) ---		1.4(0.7)/6.5 21.0 197(855) TR-20		1.6(0.8)/6.5 21.0 197(855) TR-21		1.8(1.0)/6.5 21.0 197(855) TR-62		2.0(1.1)/6.5 21.0 197(855) TR-63
196-010M 196-012M 196-014M 196-016M 196-020M 197-012M 197-014M 197-016M 197-018M 197-020M																			
圆头锥型 Round End Cone/CR/TR																			

	1.2(0.7)/8.0 22.0 198(856) ---		1.4(0.8)/8.0 22.0 198(856) TR-25		1.6(0.9)/8.0 22.0 198(856) TR-24		1.8(1.0)/8.0 22.0 198(856) TR-23		1.2(0.6)/10.0 22.0 199(850) ---		1.4(0.7)/10.0 22.0 199(850) TR-11		1.6(0.8)/10.0 22.0 199(850) TR-12		1.8(0.9)/9.0 22.0 199(850) TR-13		2.2(1.3)/9.0 22.0 199(850) TR-14		2.2(1.3)/11.0 24.0 200(850L) TR-15
198-012M 198-014M 198-016M 198-018M 199-012M 199-014M 199-016M 199-018M 199-022M 200-022M																			

	2.4(1.4)/11.0 24.0 200(850L) TR-19		1.0(0.8)/2.8 19.0 233(808) EX-41		1.4(1.1)/4.0 19.0 234(808L) EX-31		1.8(1.3)/5.0 19.0 237(808R) EX-20		2.1(1.6)/5.0 19.0 237(808R) EX-21		3.2(1.8)/5.0 19.0 237(808R) EX-26
200-024M 233-010M 234-014M 237-018M 237-021M 237-032M											
圆头倒截锥型 Round End Inverted cone/EX											

	1.2/6.0 21.0 248(861) FO-41		1.4/6.0 21.0 248(861) FO-20		1.4/8.0 22.0 249(862) FO-21		1.6/8.0 22.0 249(862) FO-22		1.2/10.0 24.0 250(863) FO-11		1.4/10.0 24.0 250(863) FO-12		1.6/10.0 24.0 250(863) FO-13		1.6/3.2 19.0 257(368,830) FO-30		1.8/3.6 19.0 257(368,830) FO-32
248-012M		248-014M		249-014M		249-016M		250-012M		250-014M		250-016M		257-016M		257-018M	
尖头圆柱型 Ogival end Cylinder/FO												花蕾芽型 Bud/FO					

	2.1/5.0 19.0 257(368,830) ---		2.3/5.0 19.0 257(368,830) FO-34		2.5/5.0 19.0 257(368,830) FO-24		2.8/5.0 19.0 257(368,830) FO-25		3.0/5.0 19.0 257(368,830) FO-26		3.2/5.0 19.0 257(368,830) FO-27		1.8/3.5 19.0 277(379,833) ---		2.0/3.9 19.0 277(379,833) TR-22		2.3/4.5 19.0 277(379,833) ---		3.0/5.0 19.0 277(379,833) TR-28
257-021M		257-023M		257-025M		257-028M		257-030M		257-032M		277-018M		277-020M		277-023M		277-030M	
蛋型 Egg Football/TR																			

	1.2/6.0 21.0 288(867,877) SO-20		1.2/8.0 22.0 289(868,878) ---		1.4/8.0 22.0 289(868,878) SO-21		1.6/8.0 22.0 289(868,878) ---		1.2/10.0 24.0 290(869,879) ---		1.4/10.0 24.0 290(869,879) ---		1.6/10.0 24.0 290(869,879) ---
288-012M		289-012M		289-014M		289-016M		290-012M		290-014M		290-016M	
鱼雷圆柱型 Torpedo/SO													

	1.6(1.0)/8.0 22.0 298(878K) ---		1.8(1.1)/8.0 22.0 298(878K) ---		1.6(1.0)/10.0 24.0 299(879K) ---		1.8(1.1)/10.0 24.0 299(879K) ---		1.4/8.0 22.0 115(842SP) BC-11		1.6/8.0 22.0 178(848SP) BC-12		1.6/5.0 19.0 465(820) ---
298-016M		298-018M		299-016M		299-018M		115-014M		178-016M		465-016M	
鱼雷圆锥型 Torpedo,Conical						邻面去釉 Interdental							

	cutting depth 0.3mm 1.6(1.0)/6.0 21.0 552(834) VP-20		cutting depth 0.4mm 1.8(1.0)/6.0 21.0 552(834) VP-21		cutting depth 0.5mm 2.1(1.1)/6.0 21.0 552(834) VP-22		 1.8(1.1)/4.5 19.0 544(845KR) RS-31 or EX-38		 1.8(1.1)/7.5 21.0 545(846KR) RS-21 or EX-29		 1.6(0.8)/8.0 22.0 546(847KR) 847KR016		 1.8(0.9)/10.0 24.0 553(848KR) 848KR018		 1.9(1.2)/10.5 23.0 553(848KR) RS-11 or EX-17
552-016M		552-018M		552-021M		544-018M		545-018M		546-016M		553-018M		553-019M	
定深车针 Depth Marking				平头圆柱型 Taper, Round Edge											


	1.6(0.9)/4.0 16.5 170(845) TF-S31		1.8(1.0)/5.0 17.0 170(845) TF-S23		2.1(1.5)/5.0 17.0 170(845) TF-S22		1.6(0.8)/6.5 18.5 171(846) TF-S21		1.4(0.7)/6.5 18.5 171(846) TF-S20		1.4(0.7)/6.5 18.5 197(855) TR-S20		1.6(0.8)/6.5 18.5 197(855) TR-S21		1.8(0.9)/9.0 21.0 199(850) TR-S13		1.4/8.0 19.5 249(862) FO-S21
170-S016M		170-S018M		170-S021M		171-S016M		171-S014M		197-S014M		197-S016M		199-S018M		249-S014M	
短柄 Short																	

	1.6(0.9)/4.0 16.5 170(845) TF-S31		1.8(1.0)/5.0 17.0 170(845) TF-S23		2.1(1.5)/5.0 17.0 170(845) TF-S22		1.6(0.8)/6.5 18.5 171(846) TF-S21		1.4(0.7)/6.5 18.5 171(846) TF-S20		1.4(0.7)/6.5 18.5 197(855) TR-S20		1.6(0.8)/6.5 18.5 197(855) TR-S21		1.8(0.9)/9.0 21.0 199(850) TR-S13		1.4/8.0 19.5 249(862) FO-S21
170-S016M		170-S018M		170-S021M		171-S016M		171-S014M		197-S014M		197-S016M		199-S018M		249-S014M	

								
Super Coarse Black	001(801) BR-41S	001(801) BR-40S	001(801) BR-31S	001(801) BR-28S	001(801) BR-26S	001(801) BR-25S	494(802L) ENDO-012S	494(802L) ENDO-014S
黑色环极粗砂	001-014S	001-016S	001-018S	001-021S	001-023S	001-025S	494-012S	494-014S
圆球型 Ball, Round /BR							根管治疗车针 /Root canal therapy	


									
494(802L) ENDO-016S	494(802L) ENDO-018S	010(805) SI-48S	038(811) EX-11S	039(811L) EX-12S	237(808R) EX-20S	237(808R) EX-21S	109(835) SF-31S	109(835) SF-42S	110(836) SF-21S
494-016S	494-018S	010-016S	038-032S	039-035S	237-018S	237-021S	109-012S	109-014S	110-014S
倒锥型 Inverted Cone /SI			双圆锥型 Double cone/EX		圆头倒截锥型 Round End Inverted cone/EX		平头圆柱型 Flat End Cylinder/SF		

									
110(836) SF-15S	111(837) SF-11S	111(837) SF-12S	111(837) SF-13S	141(881) SR-11S	141(881) SR-12S	141(881) SR-13S	165(858) TC-21S	166(859) TC-10S	166(859) TC-11S
110-016S	111-012S	111-014S	111-016S	141-012S	141-014S	141-016S	165-014S	166-014S	166-016S
圆头柱型 Round End, Cylinder/SR							尖头锥型 Pointed cone/TC		

						
166(859) TC-16S	171(846) TF-21S	172(847) TF-09S	172(847) TF-10S	173(848) TF-11S	173(848) TF-12S	173(848) TF-13S
166-018S	171-016S	172-014S	172-016S	173-014S	173-016S	173-018S
平头锥型 Flat End Cone/TF						

					
197(855) TR-21S	198(856) TR-25S	198(856) TR-24S	199(850) TR-11S	199(850) TR-12S	199(850) TR-13S
197-016S	198-014S	198-016S	199-014S	199-016S	199-018S
圆头锥型 Round End Cone/CR/TR					

						
248(861) FO-41S	248(861) FO-20S	249(862) FO-21S	249(862) FO-22S	250(863) FO-11S	250(863) FO-12S	250(863) FO-13S
248-012S	248-014S	249-014S	249-016S	250-012S	250-014S	250-016S
尖头圆柱型 Ogival end Cylinder/FO						

									
257(368,830) FO-32S	257(368,830) ---	257(368,830) FO-34S	257(368,830) FO-24S	257(368,830) FO-25S	257(368,830) FO-26S	257(368,830) FO-27S	277(379,833) TR-22S	277(379,833) TR-28S	289(868,878) SO-21S
257-018S	257-021S	257-023S	257-025S	257-028S	257-030S	257-032S	277-020S	277-030S	289-014S
花蕾芽型 Bud/FO							蛋型 Egg Football/TR		鱼雷圆柱型 Torpedo/SO

								
Coarse Green	1.2 19.0 001(801) BR-46C	1.4 19.0 001(801) BR-41C	1.6 19.0 001(801) BR-40C	1.8 19.0 001(801) BR-31C	2.1 19.0 001(801) BR-28C	2.3 19.0 001(801) BR-26C	2.5 19.0 001(801) BR-25C	2.9 19.0 001(801) BR-23C
绿色环粗砂								
圆球型 Ball, Round /BR								

								
1.2/10.0 21.0 494(802L) ENDO-012C	1.4/10.0 21.0 494(802L) ENDO-014C	1.6/10.0 21.0 494(802L) ENDO-016C	1.8/10.0 21.0 494(802L) ENDO-018C	1.4(1.0)/1.6 19.0 010(805) SI-47C	1.6(1.1)/1.8 19.0 010(805) SI-48C	3.2/5.0 19.0 038(811) EX-11C	3.5/7.0 20.0 039(811L) EX-12C	4.2/1.8 19.0 068(909) WR-13C
494-012C	494-014C	494-016C	494-018C	010-014C	010-016C	038-032C	039-035C	068-042C
根管治疗车针/Root canal therapy			倒锥型 Inverted Cone / SI		双圆锥型 Double cone/EX		轮型 Wheel/WR	

								
1.2/4.0 19.0 109(835) SF-31C	1.4/4.0 19.0 109(835) SF-42C	1.6/4.0 19.0 109(835) ---	1.2/6.0 21.0 110(836) SF-20C	1.4/6.0 21.0 110(836) SF-21C	1.6/6.0 21.0 110(836) SF-15C	1.2/8.0 22.0 111(837) SF-11C	1.4/8.0 22.0 111(837) SF-12C	1.6/8.0 22.0 111(837) SF-13C
109-012C	109-014C	109-016C	110-012C	110-014C	110-016C	111-012C	111-014C	111-016C
平头圆柱型 Flat End Cylinder/SF								

					
1.2/8.0 22.0 130(885) CP-11C	1.4/8.0 22.0 130(885) CP-12C	1.6/8.0 22.0 130(885) CP-13C	1.2/10.0 24.0 131(886) ---	1.4/10.0 24.0 131(886) ---	1.6/10.0 24.0 131(886) ---
130-012C	130-014C	130-016C	131-012C	131-014C	131-016C
Ogival end Cylinder/CP					

					
1.4/6.0 21.0 140(880) SR-21C	1.2/8.0 22.0 141(881) SR-11C	1.4/8.0 22.0 141(881) SR-12C	1.6/8.0 22.0 141(881) SR-13C	1.4/10.0 23.0 142(882) ---	1.6/8.0 23.0 142(882) ---
140-014C	141-012C	141-014C	141-016C	142-014C	142-016C
圆头柱型 Round End, Cylinder/SR					

						
1.4/8.0 22.0 165(858) TC-21C	1.4/10.0 22.0 166(859) TC-10C	1.6/10.0 22.0 166(859) TC-11C	1.8/10.0 22.0 166(859) TC-16C	1.4/11.5 24.0 167(859L) ---	1.6/11.5 24.0 167(859L) ---	1.8/11.5 24.0 167(859L) ---
165-014C	166-014C	166-016C	166-018C	167-014C	167-016C	167-018C
尖头锥型 Pointed cone/TC						

								
1.4(0.8)/4.0 19.0 170(845) TF-43C	1.4(0.7)/6.5 21.0 171(846) TF-20C	1.4(0.7)/8.0 22.0 172(847) TF-09C	1.6(0.8)/8.0 22.0 172(847) TF-10C	1.8(0.9)/8.0 22.0 172(847) TF-56C	2.3(1.3)/9.0 22.0 172(847) TF-14C	1.4(0.7)/10.0 24.0 173(848) TF-11C	1.6(0.8)/10.0 22.0 173(848) TF-12C	1.8(0.9)/10.0 22.0 173(848) TF-13C
170-014C	171-014C	172-014C	172-016C	172-018C	172-023C	173-014C	173-016C	173-018C
平头锥型 Flat End Cone/TF								

	1.4(0.7)/6.5 21.0 197(855) TR-20C		1.6(0.8)/6.5 21.0 197(855) TR-21C		1.8(1.0)/6.5 21.0 197(855) TR-62C		1.4(0.8)/8.0 22.0 198(856) TR-25C		1.6(0.9)/8.0 22.0 198(856) TR-24C		1.8(1.0)/8.0 22.0 198(856) TR-23C		1.4(0.7)/10.0 22.0 199(850) TR-11C		1.6(0.8)/10.0 22.0 199(850) TR-12C		1.8(0.9)/9.0 22.0 199(850) TR-13C		2.2(1.3)/9.0 22.0 199(850) TR-14C
197-014C 197-016C 197-018C 198-014C 198-016C 198-018C 199-014C 199-016C 199-018C 199-022C																			
圆头锥型 Round End Cone/CR/TR																			

	1.4/8.0 22.0 249(862) FO-21C		1.6/8.0 22.0 249(862) FO-22C		1.2/10.0 24.0 250(863) FO-11C		1.4/10.0 24.0 250(863) FO-12C		1.6/10.0 24.0 250(863) FO-13C		1.8/9.5 21.0 190 FO-54C
249-014C	249-016C	250-012C	250-014C	250-016C	190-018C						
尖头圆柱型 Ogival end Cylinder/FO											

	1.8/3.6 19.0 257(368,830) FO-32C		2.1/5.0 19.0 257(368,830) ---		2.3/5.0 19.0 257(368,830) FO-34C		2.5/5.0 19.0 257(368,830) FO-24C		2.8/5.0 19.0 257(368,830) FO-25C
257-018C 257-021C 257-023C 257-025C 257-028C									
花雷芽型 Bud/FO									

	1.2/8.0 22.0 289(868,878) ---		1.4/8.0 22.0 289(868,878) SO-21C		1.6/8.0 22.0 289(868,878) ---		1.2/10.0 24.0 290(869,879) ---		1.4/10.0 24.0 290(869,879) ---		1.6/10.0 24.0 290(869,879) ---
289-012C 289-014C 289-016C 290-012C 290-014C 290-016C											
鱼雷圆柱型 Torpedo/SO											

	Fine Red		1.2 19.0 001(801) BR-46F		1.4 19.0 001(801) BR-41F		1.6 19.0 001(801) BR-40F		1.8 19.0 001(801) BR-31F		2.1 19.0 001(801) BR-28F		2.3 19.0 001(801) BR-26F
红色环 较细砂 001-012F 001-014F 001-016F 001-018F 001-021F 001-023F													
圆球型 Ball, Round /BR													

	1.4/6.0 21.0 110(836) SF-21F		1.2/8.0 22.0 111(837) SF-11F		1.4/8.0 22.0 111(837) SF-12F		1.6/8.0 22.0 111(837) SF-13F		1.4/6.0 21.0 140(880) SR-21F		1.2/8.0 22.0 141(881) SR-11F		1.4/8.0 22.0 141(881) SR-12F		1.6/8.0 22.0 141(881) SR-13F		1.2/0.2 21.0 150(839) EX-18F
110-014F 111-012F 111-014F 111-016F 140-014F 141-012F 141-014F 141-016F 150-012F																	
平头圆柱型 Flat End Cylinder/SF										圆头柱型 Round End, Cylinder/SR				End Cut			

	1.2/6.0 21.0 164(852) TC-26F		1.2/8.0 22.0 165(858) ---		1.4/8.0 22.0 165(858) TC-21F		1.0/10.0 22.0 166(859) TC-27F		1.2/10.0 22.0 166(859) TC-28F		1.4/10.0 22.0 166(859) TC-10F		1.6/10.0 22.0 166(859) TC-11F		1.8/10.0 22.0 166(859) TC-16F
164-012F 165-012F 165-014F 166-010F 166-012F 166-014F 166-016F 166-018F															
尖头锥型 Pointed cone/TC															

	1.4(0.8)/4.0 19.0 170(845) TF-43F		1.4(0.7)/6.5 21.0 171(846) TF-20F		1.6(0.8)/6.5 21.0 171(846) TF-21F		1.2(0.7)/8.0 22.0 172(847) TF-08F		1.4(0.7)/8.0 22.0 172(847) TF-09F		1.6(0.8)/8.0 22.0 172(847) TF-10F		1.4(0.7)/10.0 24.0 173(848) TF-11F		1.6(0.8)/10.0 22.0 173(848) TF-12F		1.8(0.9)/10.0 22.0 173(848) TF-13F
170-014F	171-014F	171-016F	172-012F	172-014F	172-016F	173-014F	173-016F	173-018F									
平头锥型 Flat End Cone/TF																	

	1.6(0.9)/4.0 19.0 196(849) CR-12F		2.0(1.0)/5.0 19.0 196(849) CR-11F		1.4(0.7)/6.5 21.0 197(855) TR-20F		1.6(0.8)/6.5 21.0 197(855) TR-21F		1.8(1.0)/6.5 21.0 197(855) TR-62F		1.4(0.8)/8.0 22.0 198(856) TR-25F		1.6(0.9)/8.0 22.0 198(856) TR-24F		1.8(1.0)/8.0 22.0 198(856) TR-23F		1.2(0.6)/10.0 22.0 199(850) ---		1.4(0.7)/10.0 22.0 199(850) TR-11F
196-016F	196-020F	197-014F	197-016F	197-018F	198-014F	198-016F	198-018F	199-012F	199-014F										
圆头锥型 Round End Cone/CR/TR																			


	1.6(0.8)/10.0 22.0 199(850) TR-12F		1.8(0.9)/9.0 22.0 199(850) TR-13F		2.2(1.3)/9.0 22.0 199(850) TR-14F		1.4(1.1)/4.0 19.0 234(808L) EX-31F		1.8(1.3)/5.0 19.0 237(808R) EX-20F		2.1(1.6)/5.0 19.0 237(808R) EX-21F
199-016F	199-018F	199-022F	234-014F	237-018F	237-021F						
圆头倒截锥型 Round End Inverted cone/EX											







	1.2/6.0 21.0 248(861) FO-41F		1.4/6.0 21.0 248(861) FO-20F		1.4/8.0 22.0 249(862) FO-21F		1.6/8.0 22.0 249(862) FO-22F		1.2/10.0 24.0 250(863) FO-11F		1.4/10.0 24.0 250(863) FO-12F		1.6/10.0 24.0 250(863) FO-13F
248-012F	248-014F	249-014F	249-016F	250-012F	250-014F	250-016F							
尖头圆柱型 Ogival end Cylinder/FO													

	1.6/3.2 19.0 257(368,830) FO-30F		1.8/3.6 19.0 257(368,830) FO-32F		2.1/5.0 19.0 257(368,830) ---		2.3/5.0 19.0 257(368,830) FO-34F		2.5/5.0 19.0 257(368,830) FO-24F		2.8/5.0 19.0 257(368,830) FO-25F		3.0/5.0 19.0 257(368,830) FO-26F		3.2/5.0 19.0 257(368,830) FO-27F
257-016F	257-018F	257-021F	257-023F	257-025F	257-028F	257-030F	257-032F								
花蕾芽型 Bud/FO															




	1.2/6.0 21.0 288(867,877) SO-20F		1.2/8.0 22.0 289(868,878) ---		1.4/8.0 22.0 289(868,878) SO-21F		1.6/8.0 22.0 289(868,878) ---
288-012F	289-012F	289-014F	289-016F				
鱼雷圆柱型 Torpedo/SO							




	0.6 16.8 001 CD-50F		0.7/2.2 16.8 138 CD-51F		0.7/3.8 16.8 139 CD-52F		0.7/2.2 16.8 161 CD-53F		0.6/2.2 16.8 246 CD-54F		0.9/1.2 16.8 254 CD-55F		0.7/2.2 16.8 246 CD-56F		0.7/3.6 16.8 247 CD-57F		0.8/3.2 16.8 108 CD-58F		1.3/6.9 16.8 171 CD-59F
001-S006F	138-S007F	139-S007F	161-S007F	246-S006F	254-S009F	246-S007F	247-S007F	108-S008F	171-S013F										
儿童车针 FG Short Burs for Kids / CD																			





						
Extra Fine Yellow	1.2 19.0 001(801) BR-46EF	1.4 19.0 001(801) BR-41EF	1.6 19.0 001(801) BR-40EF	1.8 19.0 001(801) BR-31EF	2.1 19.0 001(801) BR-28EF	2.3 19.0 001(801) BR-26EF
黄色环 细砂	001-012EF	001-014EF	001-016EF	001-018EF	001-021EF	001-023EF
圆球型 Ball, Round /BR						






							
1.2/6.0 21.0 164(852) TC-26EF	1.2/8.0 22.0 165(858) ---	1.4/8.0 22.0 165(858) TC-21EF	1.0/10.0 22.0 166(859) TC-27EF	1.2/10.0 22.0 166(859) TC-28EF	1.4/10.0 22.0 166(859) TC-10EF	1.6/10.0 22.0 166(859) TC-11EF	1.8/10.0 22.0 166(859) TC-16EF
164-012EF	165-012EF	165-014EF	166-010EF	166-012EF	166-014EF	166-016EF	166-018EF
尖头锥型 Pointed cone/TC							

				
1.6(0.8)/6.5 21.0 171(846) TF-21EF	1.6(0.8)/8.0 22.0 172(847) TF-10EF	1.4(0.7)/10.0 24.0 173(848) TF-11EF	1.6(0.8)/10.0 22.0 173(848) TF-12EF	1.8(0.9)/10.0 22.0 173(848) TF-13EF
171-016EF	172-016EF	173-014EF	173-016EF	173-018EF
平头锥型 Flat End Cone/TF				

									
1.6(0.9)/4.0 19.0 196(849) CR-12EF	2.0(1.0)/5.0 19.0 196(849) CR-11EF	1.2(0.7)/6.5 21.0 197(855) ---	1.4(0.7)/6.5 21.0 197(855) TR-20EF	1.6(0.8)/6.5 21.0 197(855) TR-21EF	1.8(1.0)/6.5 21.0 197(855) TR-62EF	1.2(0.7)/8.0 22.0 198(856) ---	1.4(0.8)/8.0 22.0 198(856) TR-25EF	1.6(0.9)/8.0 22.0 198(856) TR-24EF	1.8(1.0)/8.0 22.0 198(856) TR-23EF
196-016EF	196-020EF	197-012EF	197-014EF	197-016EF	197-018EF	198-012EF	198-014EF	198-016EF	198-018EF
圆头锥型 Round End Cone/CR/TR									

						
1.2(0.6)/10.0 22.0 199(850) ---	1.4(0.7)/10.0 22.0 199(850) TR-11EF	1.6(0.8)/10.0 22.0 199(850) TR-12EF	1.8(0.9)/9.0 22.0 199(850) TR-13EF	1.4(1.1)/4.0 19.0 234(808L) EX-31EF	1.8(1.3)/5.0 19.0 237(808R) EX-20EF	2.1(1.6)/5.0 19.0 237(808R) EX-21EF
199-012EF	199-014EF	199-016EF	199-018EF	234-014EF	237-018EF	237-021EF
圆头倒截锥型 Round End Inverted cone/EX						

							
1.0/4.3 23.5 540(889) FO-40EF	1.2/6.0 21.0 248(861) FO-41EF	1.4/6.0 21.0 248(861) FO-20EF	1.4/8.0 22.0 249(862) FO-21EF	1.6/8.0 22.0 249(862) FO-22EF	1.2/10.0 24.0 250(863) FO-11EF	1.4/10.0 24.0 250(863) FO-12EF	1.6/10.0 24.0 250(863) FO-13EF
540-010EF	248-012EF	248-014EF	249-014EF	249-016EF	250-012EF	250-014EF	250-016EF
尖头圆柱型 Ogival end Cylinder/FO							

									
1.6/3.2 19.0 257(368,830) FO-30EF	1.8/3.6 19.0 257(368,830) FO-32EF	2.1/5.0 19.0 257(368,830) ---	2.3/5.0 19.0 257(368,830) FO-34EF	2.5/5.0 19.0 257(368,830) FO-24EF	2.8/5.0 19.0 257(368,830) FO-25EF	3.0/5.0 19.0 257(368,830) FO-26EF	3.2/5.0 19.0 257(368,830) FO-27EF	1.2/6.0 21.0 288(867,877) SO-20EF	1.4/8.0 22.0 289(868,878) SO-21EF
257-016EF	257-018EF	257-021EF	257-023EF	257-025EF	257-028EF	257-030EF	257-032EF	288-012EF	289-014EF
花蕾芽型 Bud/FO								鱼雷圆柱型 Torpedo/SO	

 25mm FG XL 316		 1.0 25.0 001(801) BR-X45	 1.2 25.0 001(801) BR-X46	 1.4 25.0 001(801) BR-X41	 1.6 25.0 001(801) BR-X40	 1.8 25.0 001(801) BR-X31	 1.2(0.9)/1.2 25.0 010(805) SI-X46	 1.2/8.0 28.0 111(837) SF-X11	 1.4/8.0 28.0 111(837) SF-X12
FG Surgical 加长柄		001-X010M	001-X012M	001-X014M	001-X016M	001-X018M	010-X012M	111-X012M	111-X014M
 1.6/8.0 28.0 111(837) SF-X13	 1.6/10.0 29.0 166(859) TC-X11	 1.4(0.7)/10.0 29.0 173(848) TF-X11	 1.6(0.8)/10.0 29.0 173(848) TF-X12	 1.8(0.9)/10.0 29.0 173(848) TF-X13	 1.6(0.8)/10.0 29.0 199(850) TR-X12	 2.2(1.3)/11.0 29.0 200(850L) TR-X15	 1.4/8.0 28.0 249(862) FO-X21	 1.6/8.0 28.0 249(862) FO-X22	
111-X016M	166-X016M	173-X014M	173-X016M	173-X018M	199-X016M	200-X022M	249-X014M	249-X016M	
 1.6 25.0 001(801) BR-X40S	 1.8 25.0 001(801) BR-X31S	 1.4 25.0 001(801) BR-X41C	 1.6 25.0 001(801) BR-X40C	 1.8 25.0 001(801) BR-X31C					
001-X016S	001-X018S	001-X014C	001-X016C	001-X018C					
 22.0mm RA, CA 204		 1.2 22.0 001(801) BR-R46	 1.4 22.0 001(801) BR-R41	 1.6 22.0 001(801) BR-R40	 1.8 22.0 001(801) BR-R31	 2.3 22.0 001(801) BR-R26	 1.0(0.8)/1.0 22.0 010(805) SI-R45	 1.2(0.9)/1.2 22.0 010(805) SI-R46	 1.4(1.0)/1.6 22.0 010(805) SI-R47
RA, CA 弯机柄		001-R012M	001-R014M	001-R016M	001-R018M	001-R023M	010-R010M	010-R012M	010-R014M
 1.6(1.1)/1.8 22.0 010(805) SI-R48	 1.0/4.0 22.0 109(835) SF-R41	 1.2/4.0 22.0 109(835) SF-R31	 1.4/4.0 22.0 109(835) SF-R42	 1.2/8.0 24.0 111(837) SF-R11	 1.4/8.0 24.0 111(837) SF-R12	 1.6/8.0 24.0 111(837) SF-R13	 1.2/8.0 24.0 141(881) SR-R11	 1.4/8.0 24.0 141(881) SR-R12	 1.4/8.0 24.0 165(858) TC-R21
010-R016M	109-R010M	109-R012M	109-R014M	111-R012M	111-R014M	111-R016M	141-R012M	141-R014M	165-R014M
 1.6/10.0 27.0 166(859) TC-R11	 1.0(0.6)/3.0 22.0 169(845) TF-R41	 1.2(0.7)/4.0 22.0 170(845) TF-R42	 1.8(1.0)/8.0 25.0 198(856) TR-R23	 1.6(0.8)/10.0 27.0 199(850) TR-R12	 1.8/3.6 22.0 257(368,830) FO-R32	 1.6 22.0 001(801) BR-R40F	 1.8 22.0 001(801) BR-R31F	 2.3 22.0 001(801) BR-R26F	 1.4/8.0 24.0 165(858) TC-R21F
166-R016M	169-R010M	170-R012M	198-R018M	199-R016M	257-R018M	001-R016F	001-R018F	001-R023F	165-R014F
 1.8/3.6 22.0 257(368,830) FO-R32F	 1.6 22.0 001(801) BR-R40EF	 1.8 22.0 001(801) BR-R31EF	 2.3 22.0 001(801) BR-R26EF	去除肉芽组织 车针套装 Degranulation Kit		 1.0 31.0	 2.5 31.0	 3.0 31.0	 3.5 31.0
257-R018F	001-R016EF	001-R018EF	001-R023EF	去除肉芽组织 车针套装 Degranulation Kit		001-R010	001-R025	001-R030	001-R035

用于制备烧结氧化锆的金刚石车针

Diamond burs for preparing sintered Zirconia

氧化锆是一种具有良好发展前景的假牙材料。氧化锆的特殊性能对制备提出了极高的要求。本系列金刚石钻头专门针对氧化锆陶瓷的极高硬度而设计。它们具有比标准金刚石钻头更长的使用寿命和更大的切削能力。它们是根管治疗中切割氧化锆冠的理想工具。

Zirconia is a material with a promising future in prosthetic dentistry, the special properties of zirconia place extreme demands on the preparation. This series of our diamond burs are tailored to the extreme hardness of zirconia ceramic. They have a long service life and much great cutting capacity than standard diamond burs. They are ideal for sectioning zirconia crowns during endodontic treatment.


									
2.1/10.0 24.0	3.2/5.0 19.0	1.2/8.0 22.0	1.4/8.0 22.0	2.3(1.3)/9.0 22.0	1.6(0.8)/10.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0
494 (802L)	038 (811)	141 (881)	141 (881)	172 (847)	173 (848)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)
Z494-021M	Z038-032M	Z141-012M	Z141-014M	Z172-023M	Z173-016M	Z290-014M	Z546-016M	Z553-018M	Z198-016M
—	ZEX-11	ZSR-11	ZSR-12	ZTF-14	ZTF-12	—	—	—	ZTR-24

									
1.4 22.0	1.2/8.0 22.0	1.4/8.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0	2.5/5.0 19.0	1.6/3.2 19.0	1.0/11.5 24.0
697 (801L)	141 (881)	141 (881)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)	257 (368,830)	277 (379,833)	167 (859L)
Z697-014F	Z141-012F	Z141-014F	Z290-014F	Z546-016F	Z553-018F	Z198-016F	Z257-025F	Z277-016F	Z167-010F
ZBR-L41F	ZSR-11F	ZSR-12F	—	—	—	ZTR-24F	ZFO-24F	—	—

									
1.4 22.0	1.2/8.0 22.0	1.4/8.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0	2.5/5.0 19.0	1.6/3.2 19.0	1.0/11.5 24.0
697 (801L)	141 (881)	141 (881)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)	257 (368,830)	277 (379,833)	167 (859L)
Z697-014EF	Z141-012EF	Z141-014EF	Z290-014EF	Z546-016EF	Z553-018EF	Z198-016EF	Z257-025EF	Z277-016EF	Z167-010EF
ZBR-L41EF	ZSR-11EF	ZSR-12EF	—	—	—	ZTR-24EF	ZFO-24EF	—	—

快速切削全陶瓷，烧结氧化锆冠的黑色金刚石车针（超粗）





Fast reducing the All-Ceramics, Zirconia crown with the Black burs (Super coarse)

									
1.2/8.0 22.0	1.4/8.0 22.0	1.6(0.8)/10.0 22.0	1.6/10.0 22.0	1.8/10.0 22.0	1.8(0.9)/9.0 22.0	2.2(1.3)/9.0 22.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	2.3/5.0 19.0
141 (881)	141 (881)	173 (848)	166 (859)	166 (859)	199 (850)	199 (850)	546 (847KR)	553 (848KR)	257 (368,830)
Z141-012S	Z141-014S	Z173-016S	Z166-016S	Z166-018S	Z199-018S	Z199-022S	Z546-016S	Z553-018S	Z257-023S
ZSR-11S	ZSR-12S	ZTF-12S	ZTC-11S	ZTC-16S	ZTR-13S	ZTR-14S	—	—	ZFO-34S

牙体制备专用车针

1. 高品质专业备牙车针，有前后牙之分。
2. 制备肩台时不会出现羽状肩台，避免了因备牙而引起的修复体无法就位和颈缘崩裂现象。
3. 保证肩台宽度一致。

Especially for All-ceramic restoration. Pin guide not Diamond coated avoids injury of the gingival edge;
508T-014XC wearing-off the contact points mesial-distal;
001-014M determination of prep-bounds
198C-016XC 021XC 508T-016XC 021XC for circular wearing-off without damaging the prep-stage;
508T-016F;021F for precise grinding of the pre-stage with the help of a precision guiding-pin.

									
1.4/8.0 21.5	1.4 19.0	1.6/8.0 20.5	1.6/9.0 21.5	1.6/9.0 21.5	1.4/8.0 21.5	1.6 16.0	2.1/8.0 19.0	2.1/9.0 20.0	2.1/9.0 20.0
508T-014XC	BR-41	198C-016XC	508T-016XC	508T-016F	508T-014XC	001-016M	198C-021XC	508T-021XC	508T-021F
508T-014XC	001-014M	198C-016XC	508T-016XC	508T-016F	508T-014XC	001-016M	198C-021XC	508T-021XC	508T-021F

Diamond Bur HP (Electroplated) 技工用金刚砂磨头 (电镀)



Sharp and Durable,
Grind or polish Ceramics
锋利、耐用
磨削平稳、不偏摆
不留黑印、适合车瓷、



Grit

- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine


Fine
Red
细砂
红圈



编号 No. FD30 FK33 FA23 FA33 FC18 FC25 FM27 FL33 FN23 FP33

This grit can grind Zirconia 可用来车氧化锆


Medium
Blue
标准砂
蓝圈



FQ23 FU40 编号 No. MD10 MD14 MD18 MD23 MD27 MD35 MD40 MK16

Round / 圆球型


Inverted Cone 倒锥型



MK21 MK27 MK37 MK47 MK60 MK80 MK100 MT60 MJ60 MA16 MA23

Wheel / 轮型


Cylinder / 平头圆柱型



MA33 MA50 MC18 MC25 MC30 MC40 MM12 MM18 MM27 MM37 MM50

Round Cylinder / 圆头柱型

Needle / 尖头锥型



MB12 MB16 MB21 MB31 MB40 ML16 ML23 ML33 ML40 ML50 MN18

Taper / 平头锥型

Round Taper / 圆头锥型

Inverted Cone 倒锥型

									
MN23	MP33	MG16	MG23	MG50	MR23	MV60	MQ23	MU50	MU100
2.3 6.0	239 3.3 8.0	250 1.6 10.0	2.3 10.0	5.0 12.0	243 2.3 5.0	266 6.0 18.0	292 2.3 14.0	304 5.0 1.5	10.0 2.5
Pearl 圓頭側錐型 Flame / 火焰					R.Taper 圓頭錐型 Torpedo 魚雷圓柱型 Wheel / 輪型				

Super Coarse Black 特粗砂黑圈									
									
001 3.5	112 5.0 11.0	141 4.0 10.0	238 7.0 12.5	251 5.0 12.0	263 7.0 12.0	266 6.0 18.0	274 6.0 14.0		
编号 No.	SD35	SA50	SC40	SP70	SG50	SV70	SV60	SF60	

Coarse Green 粗砂 绿圈											
											
	001 2.3	3.5	068 6.0 2.0	112 5.0 11.0	143 1.8 15.0	141 4.0 10.0	173 3.1 10.0	198 4.0 12.0	292 2.3 14.0	266 6.0 18.0	274 6.0 14.0
编号 No.	CD23	CD35	CJ60	CA50	CC18	CC40	CB31	CL40	CQ23	CV60	CF60

美甲用金刚砂磨头

Diamond Burs for Nails

产热少 Heatless
更锋利 Sharper
更耐磨 More Durable



型号齐 More shapes
手感好 Better feeling



USED for Nails Cutting or Finishing
More than 100 shapes for your choices
COMFORTABLE at the minimal pressure



HEATLESS no need water cooling
MORE DURABLE under our suggested PRM

红色圈 Red Ring
较细粒度 Fine

圆球型 (D)

Ball or Round
(ISO 001.514)

产品编号

Item No.



1



11D10F



11D12F



11D14F



11D16F



11D18F



11D21F



11D23F



11D25F

直径(d) x 高度(h) mm

1.0 *

1.2 *

1.4 *

1.6 *

1.8 *

2.1 *

2.3 *

2.5 *



11D27F

2.7 *



11D31F

3.1 *



11D40F

4.0 *



11D50F

5.0 *

平头圆柱型 (A1)

Cylinder
(ISO 110.514)

产品编号

Item No.



1



11A118F



11A121F



11A123F



11A125F



11A127F



11A131F

直径(d) x 高度(h) mm

1.8 x 6.0

2.1 x 6.0

2.3 x 6.0

2.5 x 6.0

2.7 x 6.0

3.1 x 6.0

圆头圆柱型 (C)

Round Cylinder
(ISO 141.514)

产品编号

Item No.



1



11C16F



11C18F



11C21F



11C23F



11C25F



11C27F



11C29F



11C31F

直径(d) x 高度(h) mm

1.6 x 8.0

1.8 x 8.0

2.1 x 8.0

2.3 x 8.0

2.5 x 8.0

2.7 x 8.0

2.9 x 8.0

3.1 x 8.0

平头锥型 (B)

Flat end
(ISO 173.514)

产品编号

Item No.



1



11B21F



11B25F



11B33F



11B40F

直径(d) x 高度(h) mm

2.1 x 10.0

2.5 x 10.0

3.3 x 10.0

4.0 x 10.0

圆头锥型 (L)

Round end Cone
(ISO 198.514)

产品编号

Item No.  1



11L16F



11L18F



11L21F



11L23F



11L31F



11L40F

直径(d) x 高度(h) mm

1.6 x 9.0

1.8 x 9.0

2.1 x 9.0

2.3 x 9.0

3.1 x 9.0

4.0 x 9.0

圆头倒截锥型 (P)

Pear
(ISO 239.514)

产品编号

Item No.  1



11P21F



11P23F



11P31F



11P40F

直径(d) x 高度(h) mm

2.1 x 7.0

2.3 x 7.0

3.1 x 7.0

4.0 x 10.0

火焰型 (R1)

Flame
(ISO 243.514)

产品编号

Item No.  1



11R116F



11R118F



11R121F



11R123F



11R125F

直径(d) x 高度(h) mm

1.6 x 8.0

1.8 x 8.0

2.1 x 8.0

2.3 x 8.0

2.5 x 8.0

长火焰型 (R2)

Long Flame
(ISO 244.514)

产品编号

Item No.  1



11R216F



11R218F



11R221F



11R223F



11R225F

直径(d) x 高度(h) mm

1.6 x 10.0

1.8 x 10.0

2.1 x 10.0

2.3 x 10.0

2.5 x 10.0

花蕾芽型 (R3)

Bud
(ISO 254.514)

产品编号

Item No.  1



11R316F



11R318F



11R321F



11R323F



11R325F

直径(d) x 高度(h) mm

1.6 x 3.0

1.8 x 4.5

2.1 x 5.0

2.3 x 5.0

2.5 x 5.0

橄榄型 (E)

Egg, Football
(ISO 277.514)

产品编号

Item No.  1



11E16F



11E18F



11E21F



11E23F

直径(d) x 高度(h) mm

1.6 x 3.6

1.8 x 4.5

2.1 x 4.5

2.3 x 4.2

蓝色圈 Bloe Ring
标准粒度 Medium

圆球型 (D)

Ball, Round
(ISO 001.524)

产品编号

Item No.  1



11D10M



11D12M



11D14M



11D16M



11D18M



11D21M



11D23M



11D25M

直径(d) x 高度(h) mm

1.0 *

1.2 *

1.4 *

1.6 *

1.8 *

2.1 *

2.3 *

2.5 *



11D27M

2.7 *



11D31M

3.1 *



11D40M

4.0 *

平头圆柱型 (A1)

Cylinder
(ISO 110.524)

产品编号

Item No.  1



11A117M



11A123M



11A125M



11A127M



11A131M

直径(d) x 高度(h) mm

1.7 x 4.5

2.3 x 6.0

2.5 x 6.0

2.7 x 6.0

3.1 x 6.0

圆头柱型 (C)

Round Cylinder
(ISO 141.524)

产品编号

Item No.  1



11C18M



11C21M



11C23M



11C25M



11C27M



11C31M

直径(d) x 高度(h) mm

1.8 x 8.0

2.1 x 8.0

2.3 x 8.0

2.5 x 8.0

2.7 x 8.0

3.1 x 8.0

平头锥型 (B)

Flat end Cone
(ISO 173.524)

产品编号

Item No.  1



11B18M



11B21M



11B23M



11B25M



11B33M



11B40M

直径(d) x 高度(h) mm

1.8 x 10.0

2.1 x 10.0

2.3 x 10.0

2.5 x 10.0

3.3 x 10.0

4.0 x 10.0

圆头锥型 (L)

Round end Cone
(ISO 198.524)

产品编号

Item No.  1



11L16M



11L18M



11L21M



11L23M



11L31M



11L40M

直径(d) x 高度(h) mm

1.6 x 9.0

1.8 x 9.0

2.1 x 9.0

2.3 x 9.0

3.1 x 9.0

4.0 x 9.0

圆头倒截锥型 (P)

Egg Footbau
(ISO 239.524)

产品编号

Item No.  1



11P21M



11P23M



11P31M



11P40M

直径(d) x 高度(h) mm

2.1 x 7.0

2.3 x 7.0

3.1 x 7.0

4.0 x 10.0

火焰型 (R1)

Flame
(ISO 243.524)

产品编号

Item No.  1



11R116M



11R118M



11R121M



11R123M



11R125M

直径(d) x 高度(h) mm

1.6 x 8.0

1.8 x 8.0

2.1 x 8.0

2.3 x 8.0

2.5 x 8.0

长火焰型 (R2)

Long Flame
(ISO 244.524)

产品编号

Item No.  1



11R216M



11R218M



11R221M



11R223M



11R225M

直径(d) x 高度(h) mm

1.6 x 10.0

1.8 x 10.0

2.1 x 10.0

2.3 x 10.0

2.5 x 10.0

花蕾芽型 (R3)

Bud
(ISO 254.524)

产品编号

Item No.  1



11R316M



11R318M



11R321M



11R323M



11R327M



11R333M

直径(d) x 高度(h) mm

1.6 x 3.0

1.8 x 4.5

2.1 x 5.0

2.3 x 5.0

2.7 x 5.0

3.3 x 6.0

橄榄型 (E)

Egg, Football
(ISO 277.524)

产品编号

Item No.  1



11E16M



11E18M



11E21M



11E23M

直径(d) x 高度(h) mm

1.6 x 4.0

1.8 x 3.6

2.1 x 4.5

2.3 x 4.2



品名	FG 十支黄色塑料盒
Name	FG 10 pcs yellow plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支红色塑料盒
Name	FG 10 pcs red plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支蓝色塑料盒
Name	FG 10 pcs blue plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支绿色塑料盒
Name	FG 10 pcs green plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支黑色塑料盒
Name	FG 10 pcs black plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 五支黄色塑料盒
Name	FG 5 pcs yellow plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支红色塑料盒
Name	FG 5 pcs red plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支蓝色塑料盒
Name	FG 5 pcs blue plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支绿色塑料盒
Name	FG 5 pcs green plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支黑色塑料盒
Name	FG 5 pcs black plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 十支推拉盒
Name	FG 10 pcs plastic box
尺寸/Size	35mm*33mm
重量/Weight	6.8g




品名	HP 30支塑料盒
Name	HP 30 pcs plastic box
尺寸/Size	144mm*71mm
重量/Weight	77.2g



品名	HP 十支塑料盒
Name	HP 10 pcs plastic box
尺寸/Size	108mm*69mm
重量/Weight	30.5g



品名	HP 十支吸塑
Name	HP 10pcs blister card
尺寸/Size	61mm*165mm
重量/Weight	22.4g



品名	FG 五支吸塑
Name	FG 5pcs blister card
尺寸/Size	95.5mm*56mm
重量/Weight	4.8g



品名	FG 五支吸塑
Name	FG 5pcs blister card
尺寸/Size	81mm*47mm
重量/Weight	4.2g




品名	HP 五支吸塑
Name	HP 5pcs blister card
尺寸/Size	98.5mm*71mm
重量/Weight	10.5g



品名	FG 十支吸塑
Name	FG 10pcs blister card
尺寸/Size	80mm*69mm
重量/Weight	6.7g



品名	FG 十支吸塑
Name	FG 10pcs blister card
尺寸/Size	105.5mm*80mm
重量/Weight	6.9g



品名	FG 五连包
Name	FG Five packs
尺寸/Size	149.5mm*50mm
重量/Weight	3.2g

产品套装 Product Kit



备冠/前牙套装
Crown Preparation Anterior Teeth Kit



备冠/后牙套装
Crown Preparation Molar Teeth Kit



备洞套装
Cavity Preparation Kit



(高) 嵌体预备套装
Inlay Preparation Kit



烤瓷牙抛光套装
Ceramics/Composite Polishing Kit



儿童微创修复套装
Kiddies Kit



烤瓷制备套装
Ceramic Restoration Kit



备肩台套装
Porcelain Shouldered Abutment Polishing Kit



烤瓷制备/修复套装
Porcelain Preparation / Repair Kit



制备/修复套装
Preparation Repair Kit



护龈综合套装
Gingival Bevel Burs Kit



前牙全冠牙体预备套装
Crown Preparation Anterior Teeth Kit



磨牙全冠牙体预备套装
Crown Preparation Molar Teeth Kit



复合树脂金刚砂车针套装
Composite Resin Kit



微创窝洞制备套装 (粗砂)
Minimally Invasive Cavity Preparation Coarse Kit



微创窝洞制备套装 (细砂)
Minimally Invasive Cavity Preparation Fine Kit



树脂修形套装
Composite Repair Kit



陶瓷复合抛光套装
Porcelain Ceramics Composite Polishing Kit



瓷冠制备套装
Porcelain Crown Preparation Kit



常规粗砂套装
Regular Coarse Sand Kit



磨牙修复套装
Molars Restoration Kit




备肩台套装
Shoulder preparation Kit




备冠/前牙套装
Crown Preparation Anterior Teeth Kit




邻面去釉套装
Interproximal enamel kit



氧化锆金刚砂车针套装
Zirconia Diamond Kit



氧化锆金刚砂车针套装
Diamond Zirconia Finishing Kit



氧化锆金刚砂车针套装
Zirconia Kit



破除氧化锆金刚砂车针套装
Diamond Zirconia Cutting Kit



去除肉芽组织车针套装
Degranulation Kit