



CAD/CAM Burs


CAD/CAM milling tools

Visit our store below:

<https://www.dentsma.com/product-category/lab-supplies/dental-milling-tools/>

Email: dentsma@163.com

Zirconia/Ti/CoCr/PMMA/Wax




CAD/CAM Burs-Zirkonzahn

Zirconia/Ti/CoCr/PMMA/Wax



Milling Equipment-Zirkonzahn

	N0.	R(mm)	RC(mm)	D(mm)	TL(mm)	Φ(mm)
Zirkonzahn-M5 Zirconia milling burs	T1	1.0	18	3	57	2.0
	T2	0.8	16	3	57	1.5
	T3	0.5	16	3	57	1.0
		0.3	10	3	57	0.6
		0.3	6	3	57	0.5
Zirkonzahn-M1 Zirconia milling burs	T1	1.0	16	6	50	2.0
	T2	0.5	16	6	50	1.0
	T3	0.3	10	6	50	0.6
		0.3	6	6	50	0.5
Zirkonzahn Mental milling burs	T1	1.5	16	6	50	3.0
	T2	1.0	16	6	50	2.0
	T3	0.5	16	6	50	1.0



CAD/CAM Burs-Amann Girrbach

Zirconia/Ti/CoCr/PMMA/Wax



	No.	R(mm)	RC(mm)	D(mm)	TL(mm)	Φ(mm)
AMANN Glass ceramic milling burs	T1	1.8	14	3	43	1.8
	T2	1.4	14	3	43	1.4
	T3	1.0	14	3	43	1.0
	T4	0.4	14	3	43	0.4
AMANN Zirconia milling burs	T1	1.25	8	3	47	2.5
	T2	0.5	16	3	47	1.0
	T3	0.3	13	3	47	0.6

TL: 总长 Overall length
RC: 有效长 flute length
D: 直径 Diameter

CAD/CAM Burs-Sirona

Zirconia/Ti/CoCr/PMMA/Wax



Milling Equipment-Sirona						
	No.	R(mm)	RC(mm)	D(mm)	TL(mm)	Φ (mm)
Sirona Zirconia milling burs	T1	1.25	20	3	48	2.5
	T2	0.5	16	3	48	1.0
	T3	0.3	10	3	48	0.6

TL:总长 Overall length

RC:有效长flute length

The diagram illustrates a dental milling bur with a cylindrical body and a tapered cutting edge. Dimension D represents the diameter of the body. Dimension TL represents the overall length from the base to the tip. Dimension RC represents the length of the flute, measured from the base to the start of the taper.

CAD/CAM Burs-VHF

Zirconia/Ti/CoCr/PMMA/Wax



	N0.	R(mm)	RC(mm)	D(mm)	TL(mm)	Φ (mm)
VHF Glass ceramic milling burs	T1	1.2	15	3	35	2.4
	T2	0.5	8	3	35	1.0
	T3	0.3	8	3	35	0.6
VHF Zirconia milling burs	T1	1.25	16	6	50	2.5
	T2	0.5	10	6	50	1.0
	T3	0.3	10	3	50	0.6

TL: 总长 Overall length
RC: 有效长 flute length

CAD/CAM Burs-IMES ICORE

Zirconia/Ti/CoCr/PMMA/Wax



Milling Equipment-IMES ICORE

Milling Equipment-IMES ICORE						
	No.	R(mm)	RC(mm)	D(mm)	TL(mm)	Φ(mm)
IMES ICORE-250i Glass ceramic milling burs	T1	1.25	15	3	42	2.5
	T2	0.5	10	3	42	1.0
	T3	0.3	10	3	42	0.6
IMES ICORE-350i Glass ceramic milling burs	T1	1.25	15	6	40	2.5
	T2	0.5	10	6	40	1.0
	T3	0.3	10	6	40	0.6
IMES ICORE-250i Zirconia ceramic milling burs	T1	1.25	15	3	42	2.5
	T2	0.5	10	3	42	1.0
	T3	0.3	10	3	42	0.6
IMES ICORE-350i Zirconia ceramic milling burs	T1	1.25	15	6	50	2.5
	T2	0.5	10	6	50	1.0
	T3	0.3	10	6	50	0.6



CAD/CAM Burs-IMES ICORE

Zirconia/Ti/CoCr/PMMA/Wax



IMES ICORE-350i
Metal milling burs

Milling Equipment-IMES ICORE

	No.	R(mm)	RC(mm)	D(mm)	TL(mm)	Φ(mm)
IMES ICORE-350i Metal milling burs	T1	1.5	15	6	50	3.0
	T2	1.0	12	6	50	2.0
	T3	0.75	10	6	50	1.5
	T4	0.5	10	6	50	1.0
	T5	1.5	14	6	50	1.5
	T6	0.5	3	6	50	0.5
	T7	1.0	5	6	50	1.0
	T8	2.0	7	6	50	2.0
	T9	2.2	18	6	50	2.2
	T10	1.5	13	6	50	1.5
	T11	2.0R0.2	5	6	50	2.0
	T12	1.5R0.1	14	6	50	1.5



CAD/CAM Burs-PM7

Zirconia/Ti/CoCr/PMMA/Wax



	No.	R(mm)	RC(mm)	D(mm)	TL(mm)	Φ (mm)
PM7 Metal milling burs	T1	1.5	15	6	50	3.0
	T2	1.0	12	6	50	2.0
	T3	0.75	10	6	50	1.5
	T4	0.3	5	6	50	0.6

TL: 总长 Overall length
RC: 有效长 flute length
RC: flute length



Zirconia/Ti/CoCr/PMMA/Wax

