

01

Anterior band removing pliers



Anterior Band Removing Pliers

Removes orthodontic bands from anterior teeth safely, without risk of enamel trauma.

02



Curved Cutting Plier

Cuts archwires intraorally, providing superior access to posterior segments.

03

End backend forming pliers



End Backend Forming Pliers

Creates precise terminal backbends distal to the buccal tube.

04

End cutting pliers



End Cutting Pliers

Cuts archwires flush to prevent posterior soft-tissue irritation.

05

Fine bending pliers
(with cutter)



Fine Bending Pliers (with Cutter)

Bends delicate archwires and trims excess wire with one instrument.

06

Fine bending pliers
(without cutter)



Fine Bending Pliers (without Cutter)

Forms precise first- and second-order bends on active archwires.

07

Fine cutter pliers



Fine Cutter Pliers

Trims fine ligature wires and elastic chains with high precision.

08

KIM



KIM Pliers

Forms uniform multi-loop edgewise archwire configurations and complex steps.

09



Straight Cutting Plier

Performs heavy-duty extraoral archwire trimming prior to initial placement.

10

Ladder pliers



Ladder Pliers

Bends repeatable, multi-tiered loops without switching instruments.

11

Ligature forming pliers



Ligature Forming Pliers

Adapts and secures stainless steel ligature wires around bracket tie-wings.


12

Posterior band removing pliers



Posterior Band Removing Pliers

Leverages tight molar bands safely with cushioned occlusal seating.



INSTRUMENT NO. 13

Torque Pliers

Grips rectangular archwires securely to apply accurate third-order torque. Completes the OrthoMaster suite as the finishing instrument for rectangular wire engagement.

Complete Instrumentation Directory

- | | |
|---|--|
| 01 Anterior Band Removing Pliers | 02 Curved Cutting Plier |
| 03 End Backbend Forming Pliers | 04 End Cutting Pliers |
| 05 Fine Bending Pliers (with Cutter) | 06 Fine Bending Pliers (without Cutter) |
| 07 Fine Cutter Pliers | 08 KIM Pliers |
| 09 Straight Cutting Plier | 10 Ladder Pliers |
| 11 Ligature Forming Pliers | 12 Posterior Band Removing Pliers |
| 13 Torque Pliers | |

Order the Complete Suite

13 precision instruments · one clinical system

ONLINE STORE

[dentsma.com/product/
orthomaster-13-piece-plier-suite](https://dentsma.com/product/orthomaster-13-piece-plier-suite)

WEBSITE

www.dentsma.com

EMAIL

dentsma@163.com