



DENTAL LABORATORY WORKSTATIONS

Professional Bench Systems for Dental Technicians

PRODUCT CATALOG

Website: <https://www.dentsma.com>

Email: dentsma@163.com

Lab Supplies · Dental Lab Equipment

Workstation Features & Specifications

Our dental laboratory workstations / benches are engineered for professional dental technicians. Every model includes the essential features required for efficient, comfortable, and precise laboratory work.

Standard Features on All Models

- **Adjustable Light Lamp**
 - **Comfortable Armrests**
 - **Auto-Retrieve Air Gun**
 - **Vacuum Device**
 - **Storage Drawers**
 - **Dustproof Lens / Shield**
 - **Storage Racks**
- Flexible, high-quality lighting positioned exactly where needed
 - Ergonomic arm supports for extended working sessions
 - Convenient air tool with automatic retraction system
 - Integrated dust extraction for a clean working environment
 - Multiple smooth-gliding drawers for tools and materials
 - Protective clear shields to keep work area clean
 - Organized storage solutions integrated into the design


Construction Materials

Work Surface: **Medium Density Fireproof Board**

Overall Structure: **Stainless Steel**

Durable, corrosion-resistant frame designed for long-term professional laboratory use.
Easy to clean and maintain in daily use.

Available Colors



Available Models

Four configurations to suit different laboratory layouts and team sizes:

1.5M
Single Person
Compact single-station

1.75M
Two Persons
Mid-size dual station

2.15M
Two Persons
Spacious dual station

Corner
Two Persons
L-shaped, 3 components

1.5M — Single Person Workstation



Blue Color Option



Grey Color Option

Technical Specifications

Type:	Single Person Workstation
Net Dimensions (L × W × H):	150 × 60 × 86 cm
Packing Dimensions:	157 × 76 × 126 cm
Net / Gross Weight:	125 kg
Work Surface:	Medium Density Fireproof Board
Frame Material:	Stainless Steel
Lighting:	Adjustable lamp included
Key Features:	Armrests, Air gun, Vacuum, Drawers, Dustproof lens

1.75M — Two Persons Workstation



Green Color Option



Red Color Option

Technical Specifications

Type:	Two Persons Workstation
Net Dimensions (L × W × H):	175 × 60 × 86 cm
Packing Dimensions:	182 × 76 × 126 cm
Net / Gross Weight:	150 kg
Work Surface:	Medium Density Fireproof Board
Frame Material:	Stainless Steel
Lighting:	Dual adjustable lamps included
Key Features:	Armrests, Air guns, Vacuum, Multiple drawers, Dustproof lenses

2.15M — Two Persons Workstation



Grey Color



Green Color



Red Color

Technical Specifications

Type:	Two Persons Workstation — Spacious
Net Dimensions (L × W × H):	215 × 60 × 86 cm
Packing Dimensions:	220 × 76 × 126 cm
Net / Gross Weight:	150 kg
Work Surface:	Medium Density Fireproof Board
Frame Material:	Stainless Steel
Lighting:	Dual adjustable lamps included
Key Features:	Armrests, Air guns, Vacuum, Multiple drawers, Dustproof lenses

Corner Design — Two Persons Workstation



L-Shaped Corner Configuration — 3 Major Components for Easy Installation

Technical Specifications

Type:	Corner / L-Shaped Two Persons Workstation
Net Dimensions:	(170 × 170) × 60 × 86 cm
Packing Dimensions:	2 × (109 × 76 × 126 cm) + 86 × 87 × 108 cm
Net / Gross Weight:	190 kg
Components:	3 major components — easy installation
Work Surface:	Medium Density Fireproof Board
Frame Material:	Stainless Steel
Lighting:	Dual adjustable lamps included
Key Features:	Armrests, Air guns, Vacuum, Drawers, Dustproof lenses, Storage racks