

Lab Case Storage Solution

Purpose-built pans, transport boxes and restoration cases for the dental laboratory bench — keeping stone casts, arch models, veneers, crowns and appliances organized, protected and traceable from prep to delivery.



Work Pan



Full Arch Transport Box



Dental Membrane Box



Veneer Storage Box



Denture Retainer Box

1 Work Pan — Large & Small, 12 Colors

2 Full Arch Model Transport Box

3 2-Inch Dental Membrane Box

4 Veneer Storage Box

5 Denture Retainer Box — 4 Sizes

CASE PAN SYSTEM

Work Pan

A detachable-lid case pan for organizing and transporting work between the operatory and the laboratory bench. An integrated document clip keeps the work order or patient routing slip attached directly to the pan, so paperwork never gets separated from the case.

FEATURES

- Available in two sizes — small and large — to match case volume
- Detachable lid with an integrated document clip for work orders and patient files
- 12 color options for color-coded case tracking across technicians or departments
- Rigid divided interior keeps models and small tools from shifting in transit

LABORATORY USAGE

- Transporting and protecting stone casts and aligner models
- Housing wax patterns and bite registrations during case build-up
- Corralling small technician tools between workstations

Paperwork travels with the case. Because the document clip is molded into the lid itself, the work order stays attached through every hand-off — prep, milling, finishing, and QC — without a separate clipboard or bag.




12 Color Options

CASE PAN SYSTEM

Work Pan — Sizes & Packing

Two case sizes cover single-unit cases up through full-mouth model sets.

CASE DIMENSIONS



CASE SIZE

Small Size

18 × 15 × 6 cm | Lid open height 21 cm



CASE SIZE

Large Size

21.5 × 19.5 × 10 cm | Lid open height 30 cm

SIZE	LENGTH	WIDTH	PAN DEPTH	HEIGHT (LID OPEN)
Small	18 cm	15 cm	6 cm	21 cm
Large	21.5 cm	19.5 cm	10 cm	30 cm


PACKING INFORMATION

SMALL SIZE		LARGE SIZE	
Packing Quantity	100 pcs / carton	Packing Quantity	50 pcs / carton
Carton Dimensions	87 × 21 × 80 cm	Carton Dimensions	64 × 49 × 68 cm
Shipping Weight	20 kg	Shipping Weight	22 kg

MODEL TRANSPORT CASE

Full Arch Model Transport Box

A rigid transport case built for full-arch stone or 3D-printed models. Snap-lock lids and high-density foam hold upper and lower arches firmly in place, while a built-in 2-inch mini case keeps small restorations traveling in the same box.



FEATURES

- Snap-lock lids keep the box securely closed in transit
- High-density foam holds upper and lower full-arch models firmly in place
- Contoured foam cutouts keep the two arches separated so models cannot strike each other
- Built-in 2-inch mini case stores small restorations, implants, or single crowns alongside the models
- Ideal for laboratories shipping finished or in-progress cases to clinics

Single package, two case types. The dedicated 2-inch compartment means full-arch models and their matching restorations or implant components can ship together in one sealed box — reducing lost parts and separate courier trips.

PACKAGING SPECIFICATIONS

Box Size (single unit)	18.5 × 8.5 × 6 cm
Carton Quantity	100 pcs per carton
Carton Size	47 × 38 × 63 cm
Gross Weight	16 kg per carton

RESTORATION PROTECTION CASES

Membrane & Veneer Storage Boxes

Compact cases designed to keep single-unit and multi-unit restorations suspended and protected during handling, polishing, and shipping.



2-Inch Dental Membrane Box

2-INCH DENTAL MEMBRANE BOX

- Two layers of highly elastic plastic film suspend single crowns, short bridges, and implants in mid-air
- The floating design keeps delicate ceramic margins from contacting the hard outer case wall during transit
- Compact 2-inch footprint fits inside larger transport boxes or shipping envelopes

Case Size 5.5 × 4.5 × 2 cm

VENEER STORAGE BOX

- Protective lip-film layers seal each compartment closed
- Non-slip compartment design keeps veneers from sliding during handling
- Precision tooth-shaped silicone slots cradle each veneer without contact damage

SIZE	DIMENSIONS
2-Inch	5.5 × 4.5 × 2 cm
4-Inch	9.8 × 5.9 × 2 cm




Veneer Storage Box — Compact & Full-Length Trays

PATIENT TAKE-HOME CASE

Denture Retainer Box

A four-size retainer and appliance case line molded in odorless PP material, with a secure snap-lock closure. Available in multiple colors, with OEM branding available on request.



FEATURES

- Four depth options sized for the full range of removable appliances
- Multi-color options for patient identification or clinic branding
- Odorless PP material, safe for daily patient handling
- Secure snap-lock closure keeps the lid sealed in a bag or pocket
- OEM service available — custom colors and printed branding on request

SIZE GUIDE

MODEL	DEPTH	DIMENSIONS	BEST FOR
Model 1 — Shallow	2.6 cm	7.7 × 6.8 × 2.6 cm	Aligners & Retainers
Model 2 — Medium	3.35 cm	7.8 × 6.85 × 3.35 cm	Mouthguards & Splints
Model 3 — Deep	4.5 cm	8.05 × 7.4 × 4.5 cm	Partials & Nightguards
Model 4 — Extra Deep	6.85 cm	9.95 × 8 × 6.85 cm	Full Dentures & Cast Models

GET IN TOUCH

One Storage System, Every Case on the Bench

From the first impression to the finished restoration, the Lab Case Storage Solution line keeps models, veneers, crowns and patient appliances organized, protected, and easy to route between technicians, departments, and clinics. Mix and match sizes and colors to match your lab's workflow.

01

Work Pan

02

Full Arch Transport
Box

03

Membrane Box

04

Veneer Storage Box

05

Denture Retainer Box

Web www.dentsma.com

Email dentsma@163.com

OEM Custom colors and branding available on request