



Medical Grade Disposable Examination Gloves Series

Premium Hygienic & Diagnostic Protection Solutions

Universal Design Features

- **Ergonomic Fit:** Ambidextrous design with straight fingers and fitted rolled cuff for tear resistance.
- **Hygienic Standard:** Non-sterile, clean, and hypoallergenic formulations to minimize skin irritation.

Product 1: Powder-Free Latex Examination Gloves



Application

Ideal for high-tactile precision tasks in dental procedures and clinical exams. Provides superior elasticity and comfort without the risk of powder contamination.

Key Parameters

Length:	28 cm $\pm$ 1.0 cm
Quality Standards:	EN 455 (CE), ASTM D6319
Material:	100% Natural Rubber Latex
Powder:	Powder-Free
Sizes:	S, M, L, XL

Product 2: Powdered Latex Examination Gloves



Application

Formulated for smooth, instant donning during fast-paced clinical routines and basic diagnostic care. Offers traditional elasticity with moisture-absorbing cornstarch powder.

Key Parameters

Length:	24 cm $\pm$ 1.0 cm
Quality Standards:	EN 455 (CE Standard)
Material:	Natural Rubber Latex
Powder:	Lightly Powdered
Sizes:	S, M, L, XL

Product 3: Nitrile Examination Gloves



Application

Engineered for latex-sensitive healthcare professionals requiring strong barrier protection against chemicals and pathogens. Delivers high puncture resistance and ergonomic tactile sensitivity.

Key Parameters

Length:	24 cm $\pm$ 1.0 cm
Quality Standards:	EN 455, EN 374
Material:	Synthetic Nitrile Rubber
Powder:	Powder-Free (Latex-Free)
Sizes:	S, M, L, XL

Comprehensive Specification Matrix

Feature	Powder-Free Latex	Powdered Latex	Nitrile
Material	Natural Latex	Natural Latex	Synthetic Nitrile
Powder Content	Powder-Free	Lightly Powdered	Powder-Free
Overall Length	28 cm ± 1.0 cm	24 cm ± 1.0 cm	24 cm ± 1.0 cm
Unit Weight	5.0 g	5.0 g	3.5 g
Sizing Options	S, M, L, XL	S, M, L, XL	S, M, L, XL
Certifications	EN 455, ASTM D6319	EN 455	EN 455, EN 374

Packaging & Logistics Specifications

Inner Box Quantity

100 pcs per box
(Pop-up interfolded)

Outer Carton Quantity

10 boxes / 1,000 pcs
per carton

Carton Dimensions

32.5 × 23 × 22 cm

Gross Carton Weight

4.3 kg

40ft HQ Capacity

4,250 cartons
per container

Contact Information

DENTSMA Dental & Medical Supplies

Website: www.dentsma.com | Email: dentsma@163.com