
CAD/CAM Materials Portfolio

High-Quality Blocks & Discs for Precision Restorations

Premium zirconia, titanium, PMMA, PEEK, wax and CEREC pin-notched blocks
for dental laboratories and chairside systems.

98 / 95 / AG Systems	CE MDR & FDA	OEM Logo Service
Full VITA Classic 16	Fast sintering options	Private label packaging

Professional · Reliable · Precise

1. Zirconia Disc



High Translucent Multilayer Zirconia Disc

- **Systems:** 98 mm, 95 mm and AG system available
- **Colors:** White, full VITA Classic 16 shades, and Bleach options
- **Multilayer:** Gradient-strength and translucency integrated premium available
- **Sintering:** Single crown & 1–3 unit short bridges — fast cycle within 2 hours; standard 17-hour cycle for large-span workpieces
- **Quality:** CE MDR approved and FDA registered

OEM / Private Label Service

Full OEM support: client logo design, custom packaging and private-label options available. We help laboratories and distributors build their own brand with consistent quality and certified materials.

Ideal for crowns, bridges and high-aesthetic restorations that demand strength and natural appearance.

2. Titanium Disc



Medical-Grade Titanium Discs

- **Diameter:** 98 mm and 95 mm available
- **Standard:** US Grade 5 (Ti-6Al-4V)

OEM / Private Label Service

Logo and packaging customization available for partners who wish to market under their own brand.

Perfect for implant abutments, frameworks and durable metal restorations requiring biocompatibility and strength.

3. PMMA Disc



Multilayer PMMA Discs for Temporary & Aesthetic Restorations

- **Multilayer:** Full VITA Classic 16 shades
- **Systems:** 98 mm, 95 mm and AG system available

OEM / Private Label Service

Custom logo and packaging options supported for private-label programs.

Excellent for temporary crowns, bridges and try-in restorations with natural color gradients.

4. PEEK Disc



High-Performance Polymer Discs

- **Colors:** Pink, Nature, White, Yellow
- **Systems:** 98 mm, 95 mm and AG system available

OEM / Private Label Service

Private logo and packaging customization available.

Ideal for frameworks, removable partial dentures and metal-free alternatives with excellent biocompatibility.

5. Wax Disc



Milling Wax Discs

- **Colors:** White and Blue
- **Systems:** 98 mm, 95 mm and AG system available

OEM / Private Label Service

Custom branding and packaging options available for partners.

Designed for precise milling of patterns used in casting and pressing techniques.

6. CEREC Pin-Notched Block



Single-Tooth Blocks for Chairside Systems

- **Materials:** Zirconia, Glass Ceramic and PMMA blocks available
- **Shades:** Full range of shade options

OEM / Private Label Service

Logo and packaging customization supported.

Compatible with CEREC and similar pin-notched systems for efficient single-unit restorations.

Why Choose Our CAD/CAM Materials

- Consistent quality and dimensional accuracy across every disc and block
- Wide system compatibility — 98 mm, 95 mm and AG systems
- Fast sintering options that improve laboratory productivity
- CE MDR approved and FDA registered zirconia materials
- Strong OEM / private-label service — logo design and custom packaging
- Reliable supply for high-volume dental laboratories and clinics

Contact us for samples, technical data sheets, OEM quotation and bulk pricing.

Professional · Reliable · Precise